

SkyHook

SkyHook is an over-the-head CPAP hose holder that enables free patient movement by keeping the hose up, out of the way and free of restrictions. The hose bracket is solid 18-8 stainless steel that ends with a 2" diameter soft "U" hose holder, allowing the CPAP hose to move freely with or without a hose cover. It can be stored flat against the wall or swivels freely and silently up to 12" from the wall, a circular acrylic disc 2" in diameter and a sturdy 7/8" thick. The clear, see-through scone blends neatly with background and attaches without tools with a transparent self-bonding adhesive to almost any surface. Just peel-off film and press in place. What's more, the bracket is easily removable and stores out of sight.

#2404 SkyHook E-Z store swivel CPAP hose bracket with see-through acrylic wall scone



Rumpled bedding and the circular grooves of a CPAP hose can lead to irritating and unwanted wake-up calls. Fortunately, now there's RingMaster, a nifty little hook-and-loop strap that neatly encircles the

CPAP hose and has a swivel clip. The swivel helps to keep the hose positioned so that it doesn't tug on your mask when you change positions during sleep. Clip attaches easily to PJs, sheets, blankets, or pillows, and

RingMaster

the loop works just as nicely on a covered hose. Comes two to the pack.

#2408-2 RingMaster CPAP hose strap
(two-pack in hang tag poly bag)

#2408-12 RingMaster CPAP hose strap
(2 per poly bag; bulk pack of 12 bags)

Hose Management Kit

Being a CPAP user just got a whole lot easier with our dandy little hose management kit consisting of our Tender Tubing CPAP hose cover that's reversible (plush fleece/smooth polyester), a set of RingMaster hose-control loops and SkyHook overhead hose-control bracket. The comfy-cozy Tender Tubing cuts tubing noise, combats condensation and is machine washable; the full-length, easy-glide zipper lets you insert the hose in a snap, or remove it just as quickly and easily to wash. And the SkyHook and RingMaster mean you can spend your time sleeping and not wrangling your CPAP hose.

- #2450 Hose Management Kit** (w/ Blue Tender Tubing)
- #2452 Hose Management Kit** (w/ Black Tender Tubing)
- #2454 Hose Management Kit** (w/ Beige Tender Tubing)



PRODUCTS NOT TO SCALE

